# Magnetic couplings, stirrer drives for chemical reactors bmd and cyclone: Absolutely tight, safe and dependable

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CUSTOMIZED SOLUTIONS FOR YOUR SPECIFIC  
PROCESS REQUIREMENTS

## Magnetic couplings Magnetkupplungen

The main difference of a magnetic stirrer drive and a conventional stirrer drive using sealing solutions such as rubber gland, lip seal or mechanical seal is the fact, that the stirrer shaft is driven without any contact to the static seal. The external and internal magnet are connected contactless through a strong magnetic field. The external magnet is located outside the pressure chamber. The internal magnet inside the pressure chamber is connected to the stirrer shaft and rotates simultaneously with the motor driven external magnet. Büchi offers magnetic couplings and magnetic stirrer drives with low to high torques, which achieve high mixing performances even in highly viscous process media.

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
## bmd magnetic coupling *bmd Magnetkupplung*

Büchi magnetic couplings hermetically seal the reactor. The drive of the internal stirrer via the magnetic coupling requires a specific motor and a gearbox, which meet customer requirements regarding viscosity of the media, stirrer speed and location/zone of installation (ATEX).

Example: Powerful stirrer drive with magnetic coupling bmd 300, ATEX version, installed on a pressure reactor versoclave®.

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## Cyclone magnetic stirrer drive *Cyclone Magnetrührantrieb*

Büchi magnetic drives incorporate a magnetic coupling with a built-in motor. The stirrer speed is controlled by its own controller unit.

Example: Powerful cyclone 300 magnetic stirrer drive with controller to control the speed of the stirrer motor.

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# cyclone – magnetic stirrer drives for lab reactors, stirred autoclaves

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Magnetic drives with 75 / 300 Ncm torque (cyclone 075 – cyclone 300)

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Büchi cyclone 075 and cyclone 300 are complete stirrer drives. The power transmission of the built-in stirrer motor takes place via an integrated magnetic coupling bmd 075 / bmd 300. The drive is screwed onto the cover plate with an M36x1.5 thread and requires no further fastening elements.

The cyclone 300 controller enables stepless stirrer speed adjustment of the motor. Stirrer speed and torque can be transferred via a RS 232 interface to bds measuring / display unit and/or to the bls software. The stirrer speed can be set electrically by a 4...20 mA signal.

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Inert magnetic drive with 250 Ncm torque: cyclone 250 i and bmd 250 i

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*cyclone 250 i, bmd 250 i*

The inert magnetic coupling bmd 250 i and magnetic stirrer drive cyclone 250 i are specially developed to avoid corrosion in processes with highly aggressive media or for processes where contact with metal is not permissible. Developed on the proven Büchi magnetic coupling design, nonmetallic materials such as ceramic and PEEK are used for all components in contact with the process media. High mechanical and chemical stability guarantee trouble-free operation even under the harshest conditions. This technology allows us to provide complete non-metallic pressure reactors with powerful high-torque and leak-proof agitators.

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Büchi AG