When performance matters

chemReactors for kilolab and pilot plant

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Customized solutions for your specific process requirements

Scale-up and production of

- Active pharmaceutical ingredients APIs
- Fine chemicals / Speciality chemicals
- Agrochemicals
- Flavors & Fragrances
- Coatings, Paints & Adhesives

Features


- Corrosion resistant
- -60°C to +200°C
- Full vacuum to 0.5 bar
- Multi-purpose reactors
- Customized solutions

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1

Agitator:

- With powerful motor and adjustable transmission or frequency converter
- Double-acting, dry-running mechanical seal
- ATEX approved for Zone 0/1

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2

Rupture disc:

- Protection of equipment and operator against overpressure

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3

Scaffolding:

- Made of stainless steel for secure installation of all components

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4

buchiflex flexible glass connection:

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Büchi chemReactors

– Multipurpose reactor systems for pilot plants and kilolabs

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Reactor systems for pilot and kilolab

Büchi chemReactors with 15 to 250 liters glass / glass lined steel reaction vessels and glass distillation overhead show reliability and best performance in the daily operation in pilot plants and kilolabs worldwide. They are tailored to different processes, temperature and pressure ranges and allow processing of most diverse chemicals.

Advantages of the Büchi reactor systems

- Swiss quality guarantees reliable and safe operations and a long operating life.
- Büchi chemReactors' smart design enables fast adaption to a wide range of chemical processes.
- Easy cleaning of the reactors allows fast product changes which result in short down times.
- Thanks to the low need for maintenance, the chemReactor is the most cost-effective solution on the market.
- High operational safety is guaranteed by compliance with international safety tests, standards and certifications.
- Modular system: feed vessel, condensers, receivers, etc. are perfectly matched for easy operation, maintenance and cleaning of the chemReactors.

Applications

- Various syntheses
- Distillation
- Phase separation in glass overhead
- Evaporation to any desired consistency
- Multi-component reactions
- Gas dispersion into liquid
- Extraction of multi-phase mixtures
- Crystallization
- Refluxing

Approvals, certificates for cGMP validations

- CE certificates of conformity (Machinery directive/PED/ATEX)
- TA-Luft certified
- FDA conform

Facts and figures

- 15–250 liter glass / glass lined steel reactors with distillation glass overhead
- -60°C (-90°C) to +200°C operating temperature
- Full vacuum to +0.5 bar overpressure

5
Impeller stirrer, baffle and bottom heating section:

- Made of glass/glass lined steel
- Corrosion resistant
- Excellent thermal process control

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6
Bottom valve:

- Flush mounted and service friendly
- Made of glass/glass lined steel

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7
Lowering device:

- Quick opening of the reactor vessel
- For easy cleaning and inspection

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8
Full visibility:

- Reactor vessel with glass dome or column for visual observation of the process

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Büchi offers three different types of chemReactors

BR, GR and CR type chemReactors differ in terms of their reactors and glass overheads.

Reactors: BR and GR reaction vessels consist of a glass reactor with a jacket for heating / cooling. CR reaction vessels consist of a glass lined bottom heating section with jacket and glass column & glass lined cover plate / glass cover and are characterized by excellent thermal conductivity, high stability and durability. The processes can be monitored visually through the reactor wall (BR/GR) or through the built-in glass column between glass lined bottom heating section and cover plate, respectively glass cover (CR). The glass lined impeller stirrer ensures excellent mixing for a wide range of processes.

Distillation glass overhead: The simple glass overhead of the BR chemReactor is designed for simple distillations in batch mode; the multi-purpose glass overhead of GR and CR chemReactors enables a variety of processes in batch and semi-batch mode.

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buchiflex – the unique flexible glass piping connection system

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Flexibility

- Safety, glass protection
- Excellent tightness (TA-Luft certified)
- Fast and easy installation
- Short maintenance measures

Flamepolished glass

- No residues
- Efficient cleaning

buchiflex system

- Self-draining piping

Metal flanges

- Safe grounding, discharge of electrostatics (ATEX)
- Flexibility through preload of the springs

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
chemReactor BR

Basic setup

- Good visibility into the reactor
- Compact distillation glass overhead

Features

- Reactor bottom: Glass
- Cover: Glass / glass lined steel cover plate
- Impeller stirrer: Glass lined steel

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chemReactor GR

Medium setup

- Good visibility into the reactor
- For a wide range of chemical processes
- Multipurpose distillation glass overhead

Features

- Reactor bottom: Glass
- Cover: Glass / glass lined steel cover plate
- Impeller stirrer: Glass lined steel

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chemReactor CR

Premium setup

- Highest return on investment (ROI)
- Largest range of chemical processes
- High mechanical stability and durability
- Optimal control of process temperature
- Multipurpose distillation glass overhead

Features

- Reactor bottom: Glass lined steel
- Reactor middle section: Glass
- Cover: Glass / glass lined steel cover plate
- Impeller stirrer: Glass lined steel

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Glass reaction vessels of Buchi chemReactors BR and GR

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BR and GR chemReactors consist of complete glass reactors with jacket (15 - 30 liters) or reactors with jacket glass vessel and glass lined cover plate (60 - 100 liters). The chemical processes can be observed through the reactor wall.

The glass lined cover plate gives the reactor the greatest possible stability - the large number of cover plate openings allows the connection of a wide range of accessories such as gas inlet tubes, pH probes, sampling systems, etc.

The reaction vessels are made of solid borosilicate glass 3.3 - for best stability and excellent heat transfer. Excellent mixing is achieved by the built-in glass lined stirrer and baffle.

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

- ① **Cover:** Glass dome / glass lined steel cover plate
- ② **Reaction vessel** in glass with heating / cooling jacket
- ③ **Bottom valve:** Dead volume free, inward opening
- ④ **Agitator:** Powerful drive with speed variation
- ⑤ **Shaft sealing:** Dry-running, double-acting
- ⑥ **Baffle** with Pt100

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Glass / glass lined steel reaction vessels of Büchi chemReactors CR

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CR reaction vessels' designs are based on reaction vessels used in chemical production. They are built for daily use in development and production and clearly stand out from laboratory reactors in terms of robustness and performance. Compared to a conventional reaction vessel made of glass lined steel, the glass top (dome, column) allows the operator to observe the processes inside the vessel. In addition, product residues are prevented from caking onto the unheated and smooth glass surface.


CR chemReactors consist of a glass lined bottom heating section with jacket and a glass dome (15 - 30 liters) or a glass lined bottom heating section with jacket and a glass column and a glass lined cover plate (60 - 250 liters). The glass lined bottom heating section provides excellent heat transfer and stable process temperatures.

The glass lined cover plate gives the reactor the greatest possible stability - the large number of cover plate openings allows the connection of a wide range of accessories such as gas inlet tubes, pH probes, sampling systems, etc. Excellent mixing for low and medium viscosity media is generated by the universal agitator system consisting of an impeller stirrer and a baffle.

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



CR162
CR252

- Borosilicate glass
Borosilikat Glas
- Glass lined steel
Emallierter Stahl

- Glass lined steel versions:**
- Standard: -60...+200°C
 - Low temperature: -90...200°C
 - Antistatic: electrically conductive

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1 Cover: Glass dome / or glass lined steel cover plate

2 Glass column

3 Glass lined, jacketed bottom heating section

4 Bottom valve: Dead volume free, inward opening

5 Agitator: Powerful drive with speed variation

6 Shaft sealing: Dry-running, double-acting

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Basic glass overhead for BR chemReactors

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Basis Setup
Basic setup

Design features

The distillation glass overheads consist of standard buchiflex components (borosilicate glass 3.3).

Büchi chemReactors are offered with three different setups:

- BR setup (basic)
- Multipurpose K (with tilted condensers)
- Multipurpose (with vertical condensers)

The basic glass overhead of the BR chemReactor is designed for simple distillations in batch mode. The multipurpose distillation glass overheads of GR and CR chemReactors allow a wide variety of processes in batch and semi-batch mode. The vacuum manifold with additional built-in valves and piping allows pressure equalization, refilling of feed vessels and emptying of receivers without affecting the process pressure / vacuum. In azeotropic distillation, the light and heavy phases can be continuously separated by using the phase separator.

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
Multipurpose glass overhead for GR/CR chemReactors

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

Options and accessories

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Vessel lift

Quick opening of the reactor allows easy cleaning (inside) and visual inspection.

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
Gas inlet tube

For gas introduction into liquid phase, nitrogen purging/blanking

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CiP spray nozzle

For effective cleaning on hard to reach areas inside closed reactors

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
pH monitoring / control

pH/Redox probe with gel electrolyte

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Conductive glass coating

Offers mechanical impact protection for glass components and complies with ATEX regulations

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Instrumentation

Anything is possible: from a simple display to a fully automated, recipe controlled system

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High containment solid feeding

High Containment Feststoff-Zugabe

Accessories for further process steps

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Mixing vessels

- 15 – 200 liter glass reactor
- Used as storage tank or mixing vessel
- Available in countless variations

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chemReactor

- 15 – 250 Liter
- -60°C bis +200°C
- Full vacuum to 0.5 bar
- buchiflex

Temperature control units from manufacturers like Huber, Lauda or Julabo can be connected to the reactors.

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Gas scrubber

- Neutralizes corrosive exhaust gases
- Suitable sizes for chemReactor, midi-Pilot and miniPilot
- Packed columns up to DN600

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3 Finish Abschluss

Nutsche filters

- Product filtration with lowering mechanism for removal of filter cake or cleaning
- 10 – 200 liter

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Customized CR60 with:


- Rectification column and reflux splitter
- pH control combined with weighing system and metering pumps
- Inert-gas (Nitrogen) blanking
- Stirrer drive with frequency converter
- Automation of different process steps
- Fully ATEX proofed

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Customer Focus

A fair partnership with our customers is paramount for us. Understanding the needs of our customers enables us to develop the best solution in each case and is the foundation of a successful cooperation for the mutual benefit.

We support our customers throughout the entire life cycle of their project, from basic design and engineering through manufacturing and installation to after-sales services and plant-optimization.

Quality and Technological Leadership

With our products and services, we strive for global quality and technology leadership. We achieve this through continuous investments in the development of the competencies of our employees, in the optimization of internal and external processes, in innovation as well as in manufacturing technologies and operating resources. Quality - Made in Switzerland.

Global Network

The maxim of our sales and services organization is «global reach – local competence». We offer advice and service in accordance with local needs and requirements through our subsidiaries and a widespread network of independent distributors.

We solely cooperate with selected distributors and put emphasis on long-standing and partnership-based relationships. The continuous training and professional development of our partners has great importance for us and forms the basis for local customer service at the highest level.

Experience and Innovation

The basis of a successful corporate development are products and services which fulfil and eventually exceed the continuously changing needs and requirements of our customers. Through ongoing innovation of our products and services, manufacturing technologies and internal value-creation processes, we deliver cutting edge technology to our customers.

Having realized hundreds of projects around the globe, we know exactly what it takes to get your reactor system built, on site, on time and performing. Use our experience to your benefit.

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