



### Setup details

Unistat® 610w & Buchi Glas Uster reactor

Temperature range: -60...200 °C  
 Cooling power:  
 7.0 kW @ 200...0 °C  
 6.4 kW @ -20 °C  
 3.3 kW @ -40 °C  
 0.8 kW @ -60 °C  
 Heating power: 6.0 kW  
 Hoses: 2x1.5 m; M38x1.5 (#6656)  
 HTF: DW-Therm (#6479)  
 Reactor: 20 litre jacketed glass pressure reactor  
 Reactor content: 15 litre DW-Therm (#6479)  
 Stirrer speed: 70 rpm  
 Control: process

## Unistat® 610w

**Cooling a Buchi Glas Uster 20-litre reactor from 100 °C to 20 °C**

### Requirement

This case study shows the temperature profile of a Unistat 610w working to cool a 20-litre reactor from 100 °C to 20 °C.

### Method

The Unistat and reactor were connected using two 1.5 m insulated metal hoses. The reactor was filled with 15 litre of "M90.055.03", a Huber supplied silicon based HTF.

### Results

The jacket temperature cools to approx. -43 °C to bring the process temperature rapidly to its set point in 40 minutes.

